

Slides, videos, links and more:

<https://github.com/physicell-training/ws2021>

# Session 3: Phenotype & Diffusion




Aneequa Sundus  @AneequaSundus

Furkan Kurtoglu  @FKurtogluSysBio

## PhysiCell Project

July 26, 2021



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project

**PhysiCell.org**

 @PhysiCell

# Agenda:

- Introduction
- Diffusion in PhysiCell(Microenvironment)
- Cell Phenotype
  - Cell Motility
  - Cell Mechanics
  - Cell Volume



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project

**PhysiCell.org**

 **@PhysiCell**

# Key parts of a PhysiCell model (1)

- **Microenvironment (stage):**
  - diffusing substrates
    - ♦ diffusion coefficient
    - ♦ decay rate
    - ♦ boundary conditions
    - ♦ Defined in XML configuration file
- **Cell Definitions (types of players):**
  - name
  - default phenotype (more on next page)
  - defined in XML configuration file

# A note about time steps

- PhysiCell is designed to account for the multiple time scales inherent to these problems, and has 3 time scales:
  - $\Delta t_{\text{diffusion}}$       diffusion, secretion, and uptake      (default: 0.01 min)
  - $\Delta t_{\text{mechanics}}$       cell movement      (default: 0.1 min)
  - $\Delta t_{\text{cell}}$       phenotype and volume changes      (default: 6 min)
- This allows some efficiency improvements: not all functions need to be evaluated at each time step.
- See the PhysiCell method paper. (Oddly, not in the User Guide (yet).)

# Microenvironment

- Boundary Conditions
  - By default, Von Neuman boundaries
  - Dirichlet's conditions and fine tuning
  - Dirichlet's nodes
- Define all substrates in the environment
  - Diffusion rate constant
  - Decay Rate

# XML config file Microenvironment

```
<microenvironment_setup>
  <variable name="oxygen" units="mmHg" ID="0">
    <physical_parameter_set>
      <diffusion_coefficient units="micron^2/min">100000.0</diffusion_coefficient>
      <decay_rate units="1/min">0.1</decay_rate>
    </physical_parameter_set>
    <initial_condition units="mmHg">38.0</initial_condition>
    <Dirichlet_boundary_condition units="mmHg" enabled="false">10.0</Dirichlet_boundary_condition>
  <Dirichlet_options>

    <boundary_value ID="xmin" enabled="true">10.0</boundary_value>
    <boundary_value ID="xmax" enabled="true">100.0</boundary_value>
    <boundary_value ID="ymin" enabled="false">100.0</boundary_value>
    <boundary_value ID="ymax" enabled="false">100.0</boundary_value>
    <boundary_value ID="zmin" enabled="false" hidden="true">100.0</boundary_value>
    <boundary_value ID="zmax" enabled="false" hidden="true">100.0</boundary_value>
  </Dirichlet_options>
```



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project

PhysiCell.org

@PhysiCell

# Microenvironment App

- <https://nanohub.org/resources/microenvnmtr>



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project

PhysiCell.org

 @PhysiCell

# Cell phenotype

- One of the most critical data elements in a PhysiCell Cell is ***phenotype***
- Hierarchically organize key behavioral elements:
  - Phenotype
    - ♦ **cycle**: advancement through a cell cycle model
    - ♦ **death**: one or more types of cell death
    - ♦ **volume**: cell's volume regulation
    - ♦ **geometry**: cell's radius and surface area
    - ♦ **mechanics**: adhesion and resistance to deformation ("repulsion")
    - ♦ **motility**: active motion (other than "passive" mechanics)
    - ♦ **secretion**: both release and uptake of chemical substrates. Interfaces with BioFVM
    - ♦ **intracellular**: used for interacting with intracellular models. See Sessions 10-11.
    - ♦ **molecular**: a place to store internalized substrates

**Documentation:** User Guide, Sec. 10



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project



[PhysiCell.org](https://PhysiCell.org)

@PhysiCell



# Phenotype: Mechanics

- **Mechanics** keeps parameters for adhesion and "repulsion"
  - Key parameter: maximum adhesion distance
    - ♦ a multiple of the cell's radius
  - (as a multiple of the cell's radius)
- Default model uses potential functions, but this can be supplemented or replaced.



**Documentation:** User Guide 11.5

# Cell definition: mechanics

```
<mechanics>
  <cell_cell_adhesion_strength units="micron/min">0.4</cell_cell_adhesion_strength>
  <cell_cell_repulsion_strength units="micron/min">10.0</cell_cell_repulsion_strength>
  <relative_maximum_adhesion_distance units="dimensionless">1.25</relative_maximum_adhesion_distance>

  <options>
    <set_relative_equilibrium_distance enabled="false"
      units="dimensionless">1.8</set_relative_equilibrium_distance>
    <set_absolute_equilibrium_distance enabled="false"
      units="micron">15.12</set_absolute_equilibrium_distance>
  </options>
</mechanics>
```

- **Options** give you some easy ways to *override* the cell-cell adhesion strength to accomplish other calibration goals:
  - **set\_relative\_equilibrium\_distance** lets you choose the equilibrium cell-cell spacing, as a multiple of the cell radius. It will automatically choose a `cell_cell_adhesion_strength` to meet your selected equilibrium spacing.
    - 2.0 would have an equilibrium spacing of 2 cell radii (radius of cell 1 + radius of cell 2). Don't exceed this!
    - 1.8 or 1.9 is more typical.
  - **set\_absolute\_equilibrium\_distance** allows you to choose this equilibrium distance in absolute (dimensional) units. This may or may not make sense as the cell changes size!

# mechanics app demo

- <https://nanohub.org/tools/trmechanics/>



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project


PhysiCell.org

 @PhysiCell

# Phenotype: Motility

$$\mathbf{v}_{\text{mot}} \sim s(b\mathbf{d}_{\text{bias}} + (1 - s)\mathbf{d}_{\text{rand}})$$

- **Motility** controls biased random migration
  - Migration speed  $s$
  - Bias direction  $\mathbf{d}_{\text{bias}}$
  - Migration bias  $0 \leq b \leq 1$ 
    - ♦ If  $b = 1$ , deterministic motion
    - ♦ If  $b = 0$ , purely Brownian motion
  - Persistence time  $T_{\text{per}}$



# Cell definition: motility

```
<motility>
  <speed units="micron/min">1</speed>
  <persistence_time units="min">1</persistence_time>
  <migration_bias units="dimensionless">.5</migration_bias>

  <options>
    <enabled>false</enabled>
    <use_2D>true</use_2D>
    <chemotaxis>
      <enabled>false</enabled>
      <substrate>substrate</substrate>
      <direction>1</direction>
    </chemotaxis>
  </options>
</motility>
```

- If "enabled" is set to false, the cell will not be motile, regardless of what speed you give it above.
- If you set use\_2D to true, then the cell restricts its motile motion to its current z-plane.
- chemotaxis allows you to use out-of-the-box chemotaxis:
  - set enabled to true to use this.
  - use "substrate" to choose which chemical factor it follows.
  - use direction = 1 to go up the gradient, and -1 to go against the gradient
  - **Important!!!** If the "substrate" does not match something defined in the microenvironment above, the initialization will fail.

# motility app demo

- <https://nanohub.org/tools/trmotility/>



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project

PhysiCell.org

 @PhysiCell

# Phenotype: Volume

- **volume** records the cell's sub-volumes:
  - nuclear and cytoplasmic
  - solid vs. fluid
  - calcified fraction
  - key parameters
- a very simple **default model** to regulate volume based on ODEs
  - Change the parameters, target values based on environment and cell state

$$\begin{aligned}\frac{dV_F}{dt} &= r_F(V_F^* - V_F) \\ \frac{dV_{NS}}{dt} &= r_N(V_{NS}^* - V_{NS}) \\ \frac{dV_{CS}}{dt} &= r_C(V_{CS}^* - V_{CS})\end{aligned}$$

**Documentation:** User Guide 11.3

# Cell definition: Phenotype: Volume

```
<volume>
  <total units="micron^3">2494</total>
  <fluid_fraction units="dimensionless">0.75</fluid_fraction>
  <nuclear units="micron^3">540</nuclear>

  <fluid_change_rate units="1/min">0.05</fluid_change_rate>
  <cytoplasmic_biomass_change_rate units="1/min">0.0045</cytoplasmic_biomass_change_rate>
  <nuclear_biomass_change_rate units="1/min">0.0055</nuclear_biomass_change_rate>

  <calcified_fraction units="dimensionless">0</calcified_fraction>
  <calcification_rate units="1/min">0</calcification_rate>

  <relative_rupture_volume units="dimensionless">2.0</relative_rupture_volume>
</volume>
```

- This gives both the steady-state "target" volume of the cell type and the initial volume for any cells you place in the simulation.
- Use the change rates to control how quickly cells move towards their target volume.
- The relative rupture volume is mostly useful to death models.
- Distinguish between State variables vs Target Parameters



# Phenotype: Volume app demo

- <https://nanohub.org/tools/volumetr>



**LUDDY**

SCHOOL OF INFORMATICS, COMPUTING, AND ENGINEERING

PhysiCell Project

PhysiCell.org

 @PhysiCell

# Next Session

- Cell Phenotype Continued
  - Cell Cycle
  - Cell Death
  - Cell Secretion and Uptake

# Funding Acknowledgements



## PhysiCell Development:

- Breast Cancer Research Foundation
- Jayne Koskinas Ted Giovanis Foundation for Health and Policy
- National Cancer Institute (U01CA232137)
- National Science Foundation (1720625)

## Training Materials:

- Administrative supplement to NCI U01CA232137 (Year 2)